

TIP GUIDE

Discover the wide range of cutting-edge ultrasonic bone surgery tips for the **NSK VarioSurg**

TOOTH EXTRACTION

Tooth extraction is possible without excessively cutting the alveolar bone around the remaining root, retaining the morphology of the alveolar bone. The extraction tips ruptures the periodontal ligament to streamline suture removal. Also effective in extracting ankylosed teeth.




Model SG17
Z305132



Model SG18R
Z305133



Model SG18L
Z305134



Root Extraction
(upper left No. 2 tooth removal)

Impacted tooth
(lower left No. 8 tooth removal)

BONE SURGERY

For chin bone harvests, notably for the mental foramen and mandibular ramus and general bone surgery.



Model SG1
Z305101



Model H-SG1
Z305151



Model SG1A
Z305138



For clinical cases that restrict movement distances around the tip.



Model SG8
Z305109



Model H-SG8
Z305155



Model SG8A
Z305139

For harvesting bone from the mandibular ramus.



Model SG14R
Z305122



Model SG14L
Z305123



Model H-SG8R
Z305156



Model H-SG8L
Z305157

For excisions in areas with inclines of projections, or which are otherwise hard to access with regular bone surgery tips.



Model SG2R
Z305102



Model SG2L
Z305103



For regions that are hard to access with regular tips, notably the posterior mandible or curved areas.



Model SG19
Z305135



For cutting curved surfaces



Model SG30
Z305137



For narrow cutting to reduce bone loss. Also suitable for corticotomies



Model SG58
Z305141



For removing tori on the posterior mandible or mandibular lingual surface

BONE SCRAPER

For flattening the alveolar crest or harvesting fragments.



Model SG3
Z305104



Model SG4
Z305105



Model SG5
Z305106



SINUS LIFT

(LATERAL WINDOW)

For bone fenestration and bone formations.



Model SG6D
Z305107



Model SG7D
Z305108



SINUS MEMBRANE DETACHMENT

For flattening the alveolar crest or harvesting fragments.



Model SG9
Z305110



Model SG10
Z305111



Model SG11
Z305112



**Buy 1 Get 1
FREE on all tips
+
FREE tip holder
when you order
3 types of tips***

*Terms and conditions apply. Cheapest tip will be supplied free of charge as part of the offer. Offer ends 31/08/2025.




mysk.co.uk nsk-shop.co.uk 0800 6341909



SOCKET LIFT


(CRESTAL APPROACH)

For crestal approach socket lifts. Enhancing the safety of procedures, being less invasive and enabling sinus lifts with short chair times. Osteotomy with ultrasonic vibration reduces the risk of damage to the Schneiderian membrane. Eliminating hammering taps through the osteotomy technique reduces patient discomfort and potential dizziness.



*Terms and conditions apply. Cheapest tip will be supplied free of charge as part of the offer.

For preparation of implant sites, adding starting points and other steps.



Model SG15A Z305124
Model SG15B Z305125



For preparation of implant sites, flattening bone surfaces and other steps.



Model SG15C Z305126
Model SG15D Z305127




Buy 1
Get 1 FREE
on all tips*

For expanding and forming implant preparation site when approaching from the alveolar crest.



Model SG16A Z305128
Model SG16B Z305129



Model SCL2D Z305182
Model SCL4D Z305184



For procedures such as osteotomy of maxillary sinus approaching from the alveolar crest and implant preparation site adjustments.



Model SCL1D Z305171
Model SCL3D Z305173
Model SCL5D Z305175



For elevating the sinus membrane with internal irrigation approaching from alveolar crest.



Model SCL1 Z305170
Model SCL3 Z305172
Model SCL5 Z305174



PERIO

For the removal of gingival biofilm.



Model P20-S Z305114



Model P25R-S Z305115
Model P25L-S Z305116



FREE tip holder
when you order
3 types of tips*


PERIO

(MAINTENANCE)

For the maintenance of implants and prostheses. Produced from special material to allow safe debridement without causing damage to implants and metals.



Model V10-S Z305117
Model V-P10 Y900184



Model V-P11R Y1002165
Model V-P11L Y1002165



Extra long double bent tips to easily access difficult to reach areas.



Model V-P12 Y1002167



A longer version of the popular V-P10 tip.

SCALING

For removal of supra-gingival tartar, tartar narrow interdental space, marginal gingival calculus.




Model G1-S Z305113

ENDODONTICS

Microscopic diamond coated tips for effective removal of dentin at the apex.



Model E30RD-S Z305118
Model E30LD-S Z305119
Model E31D-S Z305120
Model E32D-S Z305121



mysk.co.uk nsk-shop.co.uk 0800 6341909

